

DICOM Correction Proposal Form

Tracking Information - Administration Use Only	
Correction Proposal Number	CP-1527
STATUS	Final Text
Date of Last Update	2016/11/05
Person Assigned	Björn Nolte
Submitter Name	Renate Höcker Renate.Hoecker@siemens.com
Submission date	2015/04/09

Correction Number	CP-1527
Log Summary: Correction of X-Ray Frame Acquisition Sequence Attributes Description	
Type of Modification Addition	Name of Standard PS 3.3 - 2016d
<p>Rationale for Correction:</p> <p>C.8.19.6.8 X-Ray Frame Acquisition Macro is referenced in</p> <ul style="list-style-type: none"> a) Table A.47-2. Enhanced XA Image Functional Group Macros, b) Table A.48.4 Enhanced XRF Image Functional Group Macros c) Table A.74-2. Breast Projection X-Ray Image Functional Group Macros <p>But the description in C.8.19.6.8.1 X-Ray Frame Acquisition Sequence Attributes does not fit to c) Breast Projection X-Ray Image.</p>	
<p>Sections of documents affected</p> <p>PS 3.3, C.8.19.6.8.1</p>	
<p>Correction Wording:</p> <p>Add reference to the module referenced in Breast Projection X-Ray Image containing KVP and X-Ray Tube current in mA</p>	
<p><i>Update PS 3.3, C. C.8.19.6.8.1 as follows:</i></p>	

C.8.19.6.8.1 X-Ray Frame Acquisition Sequence Attributes

These attributes may only be used if the information is available on a frame-by-frame base. The average values for these attributes of all frames shall be stored in the same attribute in the XA/XRF Acquisition Module (Section C.8.19.3) **or in the Enhanced Mammography Image Module (Section C.8.31.1).**