## **DICOM Correction Proposal Form**

May 18, 2001

Status: Assigned

Tracking Information - Administration Use Only			
Correction Proposal Number	CP-xxx		
STATUS	Proposed		
Date of Last Update	2001/05/18		
Person Assigned	Andrei Leontiev andrei_leontiev@idx.com		
Submitter Name	Andreas Thumser Andreas.Thumser@med.siemens.de Niki Wirsz Nikolaus.Wirsz@med.siemens.de		
Submission date	2000/06/13		

Correction Number CF	CP-207	
Log Summary: Retire Unused DICOM Attributes		
Type of Modification	Name of Standard	
Correction	PS 3.6 - 2000	

## Rationale for Correction:

There are a couple of DICOM Attributes which can be found only in Part 6 (DICOM Dictionary) of the Standard:

Therapy Type (0018,0037) – introduced in Supplement 4, presumably for Therapy Module Mask Pointer(s) (0028,6030) – introduced in Supplement 4, presumably for Frame Pointers Module

Distance Source to Support (0040,0307) – introduced in Supplement 32, presumably for X-Ray Acquisition Dose Module

Referenced Procedure Step Sequence (0040,0330) – introduced in Supplement 17, intent not clear

These attributes are not referenced in any other parts of the DICOM Standard.

Sections of documents affected

PS 3.6 - 2000

Section 6 Registry of DICOM Data Elements

**Correction Wording:** 

## Remove the following entries in PS 3.6 - 2000, Section 6 Registry of DICOM Data Elements

Tag	Name	VR	VM
(0018,0037)	Therapy Type	<del>CS</del>	4
<del>(0028,6030)</del>	Mask Pointer(s)	<del>US</del>	<del>1-n</del>
(0040,0307)	Distance Source to Support	ÐS	4
(0040,0330)	Referenced Procedure Step Sequence	<del>SQ</del>	4