

DICOM Correction Item

| | | |
|---|------------------|--------|
| Correction Number | | CP-254 |
| Log Summary: Icon Image Sequence in DICOMDIR for RT Dose and Stored Print | | |
| Type of Modification | Name of Standard | |
| Addition | PS 3.3 2000 | |
| Rationale for Correction | | |
| RT Dose and Stored Print directory records contain Icon Image Sequences in DICOMDIR but are not mentioned in F.7. | | |
| Sections of documents affected | | |
| PS 3.3 F.7 | | |
| Correction Wording: | | |

F.7 ICON IMAGE KEY DEFINITION

An Icon Image may be used as a key representative of an Image, **RT Dose, Stored Print, Series or Overlay** in a corresponding Directory Record to **allow an application to display icons allowing the user to** select one or more **Image, Series or Overlay from** amongst several of them. It is based on the general purpose Image Pixel Module (See Annex C).