

DICOM Correction Item

Correction Number		CP-637
Log Summary: Additions to Anatomical Terms and Projections		
Type of Modification	Name of Standard	
Clarification	PS 3.16 2006	
Rationale for Correction		
<p>Based on the 10th Edition of Merrill's Atlas of Radiographic Poistions & Radiologic Procedures, there are several Anatomical Regions and Eponymous Names that are not included in PS3.16 2004 that should be in the standard, plus some additional corrections.</p> <p>Two typos are found in CID 4010.</p> <ul style="list-style-type: none">• acanthoparietal should be spelled acanthioparietal• parietoacanthal should be spelled parietoacanthial <p>Add views to CID 4010.</p> <ul style="list-style-type: none">• "latero-medial oblique" is already in the standard, but not in table CID 4010.• "medio-lateral oblique" is already in the standard, but not in table CID 4010. <p>For DX View Modifier in CID4011:</p> <ul style="list-style-type: none">• "lordotic" is listed in CID 4011 as a DX View Modifier. This is a patient postion and should be moved to CID 20.• "Crosstable" needs to be added as a view modifier• "Closed-Mouth" is requested to be added as a modifier. <p>For Anatomic region in CID4031:</p> <ul style="list-style-type: none">• "Eye" is already in the standard, but not in table CID 4031.• "Ilium" is needed since it is a pair body part like the Hip which is part of the pelvis. Due to the paired structure, pelvis is not an adequate description. <p>Add the name views listed in the correction section.</p>		
Sections of documents affected		
PS 3.16		
Correction Wording:		

Note: Concepts with no code value have been submitted to SNOMED, and the codes will be assigned prior to final text

Add new term to CID 20:

CID 20 Patient Orientation Modifier

**Context ID 20
 Patient Orientation Modifier**

Type: Extensible Version: 20020904-2005mmdd

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
...		
<u>SRT</u>	<u>R-40799</u>	<u>lordotic</u>

Correct Spelling and add terms in CID 4010:

CID 4010 DX View

**Context ID 4010
 DX View**

Type: Extensible Version: 20040322-2005mmdd

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
...		
SNM3	R-102C5	parietoacanthal <u>parietoacanthial</u>
SNM3	R-102C6	acanthoparietal <u>acanthioparietal</u>
...		
<u>SRT</u>	<u>R-10230</u>	<u>latero-medial oblique</u>
<u>SRT</u>	<u>R-10226</u>	<u>medio-lateral oblique</u>

Delete entry in CID 4011:

CID 4011 DX View Modifier

**Context ID 4011
 DX View Modifier**

Type: Extensible Version: 20040322-2005mmdd

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
...		
<u>SNM3</u>	<u>R-40799</u>	<u>Lordotic</u>

Add entry in CID 4011:

CID 4011 DX View Modifier

**Context ID 4011
DX View Modifier**

Type: Extensible Version: ~~20040322~~ 2005mmdd

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
...		
<u>SRT</u>		<u>Crosstable</u>

Add new terms to CID 4031:

CID 4031 Common Anatomic Regions

**Context ID 4031
Common Anatomic Regions**

Type: Extensible Version: ~~20020904~~ 2005mmdd

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
...		
<u>SRT</u>	<u>T-AA000</u>	<u>Eye</u>
<u>SRT</u>	<u>T-12340</u>	<u>Ilium</u>

Add new terms to CID 4012:

CID 4012 Projection Eponymous Name

Context ID 4012
Projection Eponymous Name
Type: Extensible Version: ~~20040322~~ 2005mmdd

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
...		
<u>SRT</u>		<u>Apple</u>
<u>SRT</u>		<u>Burman</u>
<u>SRT</u>		<u>Colcher-Sussman</u>
<u>SRT</u>		<u>Danelius-Miller</u>
<u>SRT</u>		<u>Fisk</u>
<u>SRT</u>		<u>Folio</u>
<u>SRT</u>		<u>Garth</u>
<u>SRT</u>		<u>Eraso Modification</u>
<u>SRT</u>		<u>Hirtz Modification</u>
<u>SRT</u>		<u>Kite</u>
<u>SRT</u>		<u>Moore</u>
<u>SRT</u>		<u>Neer</u>
<u>SRT</u>		<u>Rafert</u>
<u>SRT</u>		<u>Rafert-Long</u>
<u>SRT</u>		<u>Robert</u>
<u>SRT</u>		<u>Rosenberg</u>
<u>SRT</u>		<u>Stryker</u>
<u>SRT</u>		<u>Wolf</u>