

DICOM Correction Proposal

| | |
|---------------------|--|
| STATUS | New |
| Date of Last Update | 2013/06/11 |
| Person Assigned | Andrei Leontiev (andrei.leontiev@ge.com) |
| Submitter Name | IHE Cardiology, Harry Solomon |
| Submission Date | 2013/08/06 |

| | |
|---|---------|
| Correction Number | CP-1320 |
| Log Summary: Update to CID 3411 Intracoronary Devices | |
| Name of Standard PS 3.16 2011 | |
| Rationale for Correction CID 3411 was derived from a set of concepts in the National Cardiovascular Data Registry™ 2.0. IHE Cardiology is now developing a profile based on NCDR 4.x, and this context group needs to be updated with the current classifications, the addition of SNOMED Concept ID reference codes, and correction of the concept “Laser” (device, rather than light). Due to the proliferation of devices and trade names, the trade name column is deleted. <i>The final content of this CP is in discussion with NCDR and ACC Cath Report Committee, and will be finalized prior to ballot.</i> | |
| Correction Wording: | |

CID 3411 Intracoronary Devices

Context ID 3411 Intracoronary Devices

Type: Extensible

Version: 20110609 yyymmdd

| Coding Scheme Designator | Code Value | Code Meaning | Equivalent SNOMED-CT Concept ID | NCDR Equivalent | Trade Name (Informative) |
|--------------------------|-----------------------|---|---------------------------------|------------------|--------------------------|
| SRT | A-26912 | Percutaneous Transluminal Angioplasty Balloon | <u>102319006</u> | 113-1 | |
| SRT | R-002F0 | Cutting Balloon Angioplasty (CBA) Device | <u>371794009</u> | 113-2 | |
| SRT | A-25500 | Stent | <u>65818007</u> | 113-3 | |
| SRT | <u>A-25502</u> | <u>Metal Stent</u> | <u>257363003</u> | | |
| SRT | <u>A-040F8</u> | <u>Drug Coated Stent</u> | <u>411114003</u> | | |
| SRT | <u>A-040F9</u> | <u>Drug Eluting Stent</u> | <u>411191007</u> | | |
| | | <u>Covered Stent</u> | | | |
| SRT | <u>A-25600</u> | <u>Atherectomy Device</u> | <u>102312002</u> | | |
| SRT | R-002FD | Directional Coronary Atherectomy (DCA) Device | <u>371796006</u> | 113-4 | |
| SRT | A-25610 | Rotational Atherectomy Device | <u>102313007</u> | 113-5 | Rotablator™ |

| | | | | | |
|------------|--------------------------------------|---|------------------|---------------|------------------|
| SRT | R-0036F | Saline Thrombectomy Device | <u>371797002</u> | 113-6 | AngioJet™ |
| SRT | A-26920 | Transluminal Extraction Catheter (TEC) | <u>21870002</u> | 113-7 | |
| SRT | A-81080 <u>A-23000</u> | Laser Device | <u>122456005</u> | 113-8 | |
| SRT | A-10002 | Diagnostic Device | 394866000 | | |
| SRT | R-00312 | Intravascular Ultrasound (IVUS) Device | <u>371795005</u> | 113-9 | |
| SRT | R-00310 | Intracoronary Doppler guide wire | <u>371788001</u> | 113-10 | Flowire™ |
| SRT | R-00311 | Intracoronary pressure guide wire | <u>371789009</u> | 113-11 | |
| SRT | A-00D87 | Intravascular Optical Coherence Tomography Device | <u>445282004</u> | | |
| SRT | A-040ED | Brachytherapy Device | <u>228748004</u> | | |
| SRT | R-00361 | Radiofrequency Ablation Device | <u>371791001</u> | | |
| SRT | A-00927 | Guide Wire | <u>272224001</u> | | |
| SRT | A-26802 | Guiding Catheter | <u>102317008</u> | | |
| | | Chronic Total Occlusion (CTO) Catheter | | | |
| SRT | A-14610 | Vascular filter | 48295009 | | |