

1	Status	Final Text
2	Date of Last Update	2018/06/05
3	Person Assigned	David Clunie
4		mailto:dclunie@dclunie.com
5	Submitter Name	Regis Deshayes
6		mailto:regis.deshayes@zeiss.com
7	Submission Date	2016/12/07

8	Correction Number CP-1761	
9	Log Summary: Correct Acquisition Context and DateTime in Confidentiality Profile Attributes	
10	Name of Standard	
11	PS3.15 2018b	
12	Rationale for Correction:	
13	There are are two inconsistencies in Part 15's Basic Application Level Confidentiality Profile:	
14	<ul style="list-style-type: none"> • The Action Code for the Acquisition Context Sequence (0040,0555) is X whereas this is a Type 2 sequence in Part3's Table C.7.6.14-1. Acquisition Context Module Attributes, which is a mandatory Module in many IODs; its Action Code should be X/Z 	
15		
16	<ul style="list-style-type: none"> • The Action Code for the Acquisition DateTime attribute (0008,002A) is X/D whereas this is a Type 2 attribute in Part3's Table C.24-2. Encapsulated Document Module Attributes; its Action Code should be X/Z/D 	
17		
18	Correction Wording:	

Amend DICOM PS3.15 as follows (changes to existing text are bold and underlined for additions and ~~struckthrough~~ for removals):

Table E.1-1. Application Level Confidentiality Profile Attributes

Attribute Name	Tag	Retired (from ????)	In Std. Comp. IOD (from ????)	Basic Profile	Retain Safe Private Option	Retain UIDs Option	Retain Device Ident. Option	Retain Patient Chars. Option	Retain Long. Full Dates Option	Retain Long. Modif. Dates Option	Clean Desc. Option	Clean Struct. Cont. Option	Clean Graph. Option
Acquisition Context Sequence	(0040,0555)	N	Y	<u>X/Z</u>								C	
Acquisition DateTime	(0008,002A)	N	Y	<u>X/Z/D</u>					K	C			