

DICOM Correction Proposal

Status	Final Text
Date of Last Update	2023/01/21
Person Assigned	Christof Schadt
Submitter Name	Ulrich Busch (ulrich.busch@bluewin.ch)
Submission Date	2022/03/31

Correction Number	CP-2246
Log Summary:	Correct Radiotherapy Common Instance Specification
Name of Standard	PS3.3 2022e
Rationale for Correction:	<p>In C.36.4 Radiotherapy Common Instance Module the Items of Author Identification Sequence (3010,0019) invoke the Identified Person or Device Macro, Table C.17-3b. At invocation the attribute description says: "Context Group may be defined in the Macro invocation". However, this provision is useful only in a Macro itself. Here the Macro is included in a Module, and therefore the location to declare the Context Group is the specification of IODs containing this Module. The correction of that statement will help readers to find the specification referred to, which is usually located in Part 3, Appendix A in such cases.</p> <p>The text is changed accordingly. Further on the Macro title is added.</p>
Correction Wording:	

In PS3.3, modify Section C.36.4 Radiotherapy Common Instance Module

C.36.4 Radiotherapy Common Instance Module

Table C.36.4-1 specifies the Attributes that identify and describe general information in RT Second Generation IODs.

Table C.36.4-1. Radiotherapy Common Instance Module Attributes

Attribute Name	Tag	Type	Attribute Description
Instance Creation Date	(0008,0012)	1	Date the SOP Instance was created.
Instance Creation Time	(0008,0013)	1	Time the SOP Instance was created.
Content Date	(0008,0023)	1	The date the content creation started.
Content Time	(0008,0033)	1	The time the content creation started.
Author Identification Sequence	(3010,0019)	2	<p>The person that created the entire clinical content of this document.</p> <p>Zero or more Items shall be included in this Sequence.</p>
<p>>Include Table C.17-3b "<u>Identified Person or Device Macro Attributes</u>"</p>			<p>Context Group may be defined in the Macro invocation by the IOD including this Module.</p>

Attribute Name	Tag	Type	Attribute Description
			<i>The Observer Type (0040,A084) shall be PSN.</i>
Instance-Level Referenced Performed Procedure Step Sequence	(3010,0044)	1C	<p>Uniquely identifies the Performed Procedure Step SOP Instance that resulted in creation of this Instance (e.g. a Modality or Unified Procedure Step SOP Instance).</p> <p>Required if this Instance has been created as a result of a procedure step request and the Referenced Performed Procedure Step Sequence (0008,1111) is not present.</p> <p>One or more Items shall be included in this Sequence.</p>
<i>>Include Table 10-11 "SOP Instance Reference Macro Attributes"</i>			