

DICOM Correction Item

Correction Number		CP-254
Log Summary: Icon Image Sequence in DICOMDIR for RT Dose and Stored Print		
Type of Modification	Name of Standard	
Addition	PS 3.3 2000	
Rationale for Correction		
RT Dose and Stored Print directory records contain Icon Image Sequences in DICOMDIR but are not mentioned in F.7.		
Sections of documents affected		
PS 3.3 F.7		
Correction Wording:		

F.7 ICON IMAGE KEY DEFINITION

An Icon Image may be used as a key representative of an Image, **RT Dose, Stored Print, Series or Overlay** in a corresponding Directory Record to **allow an application to display icons which enable a user to** select one or more **Image, Series or Overlay from** amongst several of them. It is based on the general purpose Image Pixel Module (See Annex C).