

DICOM Correction Item

Correction Number CP-256	
Log Summary: Add new view terms to NM View Code sequence	
Type of Modification Addition	Name of Standard PS 3.16 2001
Rationale for Correction Three SNOMED codes have been identified for use in describing cardiac specific image views. These need to be added to the list of terms that can be used in the NM View Code Sequence, i.e. CID 26	
Sections of documents affected In PS 3.16, CID 26	
Correction Wording:	

Add the following 3 terms to CID 26:

Context ID 26

Nuclear Medicine Projections

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SNM3	G-A186	Short Axis
SNM3	G-A18A	Vertical Long Axis
SNM3	G-A18B	Horizontal Long Axis