

DICOM Correction Item

Correction Number CP-360	
Log Summary: Add higher precision X-ray Dose attributes	
Type of Modification Addition	Name of Standard PS 3.3 2003
Rationale for Correction Add higher precision X-ray Dose attributes (used in DX) to existing XA IODs	
Sections of documents affected PS 3.3 Section C.8.7.2	
Correction Wording:	

Table C.8-27
X-RAY ACQUISITION MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description
...			
X-Ray Tube Current in μA	(0018,8151)	3	X-Ray Tube Current in μA.
Exposure Time in μS	(0018,8150)	3	Duration of X-Ray exposure in μsec.