

## DICOM Correction Proposal Form

Correction Number		CP-365
Log Summary: Laterality in Enhanced MR		
Type of Modification	Name of Standard	
Clarification	PS 3.3 2003	
Rationale for Correction		
<p>The function of Laterality in the General Series Module is replaced in the Enhanced MR object by Frame Laterality in the Frame Anatomy functional group macro. Correct the condition on Laterality accordingly, in the same manner as was done for Image Laterality for the DX family of objects.</p>		
Sections of documents affected		
PS 3.3 C.7.3.1		
Correction Wording:		

**Table C.7-5a**  
**GENERAL SERIES MODULE ATTRIBUTES**

Attribute Name	Tag	Type	Attribute Description
...	...	...	...
Laterality	(0020,0060)	2C	<p>Laterality of (paired) body part examined. Required if the body part examined is a paired structure and Image Laterality (0020,0062) <b>or Frame Laterality (0020,9072) are is</b> not sent. Enumerated Values:</p> <p style="padding-left: 40px;">R = right L = left</p> <p>Note: Some IODs support Image Laterality (0020,0062) at the Image level <b>or Frame Laterality(0020,9072) at the Frame level in the Frame Anatomy functional group macro</b>, which can provide a more comprehensive mechanism for specifying the laterality of the body part(s) being examined.</p>