

### DICOM Correction Item

Correction Number	CP-391
Log Summary: Corrections to Mammography CAD Context Groups	
Type of Modification	Name of Standard
Correction	PS 3.16 2003
Rationale for Correction	
The name of Context ID 6018 does not correctly reflect the entire list of terms, and the code meaning values in 6019 are not all correct. The code meaning values were not assigned correctly when these terms were added to SNOMED.	
Sections of documents affected	
PS 3.16, Annex A, B	
Correction Wording:	
Update the following tables in PS 3.16 Annex A, B:	

#### TID 4005 Mammography CAD Composite Feature Body Template

The details of a composite feature are expressed in this template. It is applied to Mammography CAD Composite Feature (TID 4004).

**TID 4005**  
**MAMMOGRAPHY CAD COMPOSITE FEATURE BODY**  
**Type: Non-Extensible**

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
17			CODE	(111014, DCM, "Clockface or region")	1	U		DCID (6018) "Clockface Location <u>or</u> Region"

#### CID 6018 Clockface Location or Region

**Context ID 6018**  
**Clockface Location or Region**

**Type: Extensible Version: 20020904**

Note: In future extensions, Clockface Location or Region terms that are not derived from BI-RADS® should be added to this context group.

Coding Scheme Designator (0008,0102)	Coding Scheme Version (0008,0103)	Code Value (0008,0100)	Code Meaning (0008,0104)
			Include CONTEXT GROUP 6019

**CID 6019 Clockface Location or Region from BI-RADS®**

**Context ID 6019  
Clockface Location or Region from BI-RADS®**

**Type: Extensible Version: 20020904**

Note: From BI-RADS® 3.1, with Addendum 3.1 (National Mammography Database, E96)

<b>Coding Scheme Designator (0008,0102)</b>	<b>Coding Scheme Version (0008,0103)</b>	<b>Code Value (0008,0100)</b>	<b>Code Meaning (0008,0104)</b>
SRT	<del>1.1</del>	F-01781	1 o'clock position
SRT	<del>1.1</del>	F-01782	2 o'clock position
SRT	<del>1.1</del>	F-01783	3 o'clock position
SRT	<del>1.1</del>	F-01784	4 o'clock position
SRT	<del>1.1</del>	F-01785	5 o'clock position
SRT	<del>1.1</del>	F-01786	6 o'clock position
SRT	<del>1.1</del>	F-01787	7 o'clock position
SRT	<del>1.1</del>	F-01788	8 o'clock position
SRT	<del>1.1</del>	F-01789	9 o'clock position
SRT	<del>1.1</del>	F-0178A	10 o'clock position
SRT	<del>1.1</del>	F-0178B	11 o'clock position
SRT	<del>1.1</del>	F-0178C	12 o'clock position
SRT	<del>1.1</del>	F-0178D	Subareolar <b>positionregion</b>
SRT	<del>1.1</del>	F-0178E	Axillary tail <b>positionregion</b>
SRT	<del>1.1</del>	F-0178F	Central <b>portionregion</b> of breast <b>position</b>
SRT	<del>1.1</del>	F-01794	Axilla <b>positionregion</b>