

### DICOM Correction Item

Correction Number		CP-444
Log Summary: Clarify the use of Reference Series in Presentation State		
Type of Modification	Name of Standard	
Clarification	PS3.3 2003	
Rationale for Correction		
<p>There is no specification in the description of the Referenced Series Sequence (0008,1115) attribute to which series and instances the reference is made. Especially with the Series Instance UID there is confusion what values has to be supplied (experience of IHE Connectathon).</p> <p>This CP adds text to the relevant attribute descriptions to prevent wrong usage.</p> <p>The redundant Type 1C conditions are removed (CP-312).</p> <p>The Retrieval attributes are not considered of any value and will be removed.</p>		
Sections of documents affected		
PS 3.3 Section C.11.10		
Correction Wording:		

**C.11.10 Presentation State Module**

Table C.11.10-1 contains Attributes that identify and describe a Presentation State.

**Table C.11.10-1  
 PRESENTATION STATE MODULE ATTRIBUTES**

Attribute Name	Tag	Type	Attribute Description
...			
Referenced Series Sequence	(0008,1115)	1	Sequence of <b>Repeating</b> Items where each Item includes the Attributes of one <b>or more</b> Series <b>to which the Presentation applies</b> . <b>One or more Items shall be present.</b>
>Series Instance UID	(0020,000E)	1C	Unique identifier of a Series that is part of <b>this the Study to which the Presentation applies, defined by Study Instance UID (0020,000D)</b> . <b>Required if sequence item is present.</b>
>Retrieve AE Title	(0008,0054)	3	<b>Title of the DICOM Application Entity where the Image(s) may be retrieved on the network.</b>
>Storage Media File-Set ID	(0088,0130)	3	<b>The user or implementation specific human readable identifier that identifies the Storage Media on which the Image(s) reside.</b>
>Storage Media File-Set UID	(0088,0140)	3	<b>Uniquely identifies the Storage Media on which the Image(s) reside.</b>
>Referenced Image Sequence	(0008,1140)	1C	Sequence of <b>Repeating</b> Items where each Item provides reference to a selected set of Image SOP Class/SOP Instance pairs <b>to which the Presentation applies</b> that are part of <b>this the Study defined by Study Instance UID (0020,000D)</b> and the Series defined by Series Instance UID (0020,000E). <b>Required if a sequence item is present.</b>
>>Referenced SOP Class UID	(0008,1150)	1C	Uniquely identifies the referenced SOP Class. <b>Required if sequence item is present.</b> Shall be the same for all Images referenced by this presentation state.
>>Referenced SOP Instance UID	(0008,1155)	1C	Uniquely identifies the referenced SOP Instance. <b>Required if sequence item is present.</b>
>>Referenced Frame Number	(0008,1160)	1C	Identifies the frame numbers within the referenced SOP Instance to which the presentation applies. Note: This Attribute may be multi-valued. Required if <b>sequence item is present and</b> the referenced SOP Instance is a multi-

			frame image and the presentation does not apply to all frames.
...			