

DICOM Correction Item

Correction Number		CP-512
Log Summary: New SOP Classes added to PS 3.2 SOP Class Table		
Type of Modification	Name of Standard	
Addition	PS 3.2 - 2004	
<p>Rationale for Correction</p> <p>PS 3.2 Table A.1-2 provides guidelines how a SOP Class belongs to a category of network services (Table A.1-1).</p> <p>New SOP Classes introduced since Part P3.2 was renewed by Supplement 64 are not included in this table. This CP is adding those new SOP Classes.</p> <p>For new and supplements in progress that add new SOP Classes the update of this table should be included in the supplement.</p> <p>Also, in Table A.1-2 the SOP Class UID 1.2.840.10008.5.1.1.17.376 is included. This is not a SOP Class, but a Well-known Printer SOP Class. This entry must be removed.</p>		
Sections of documents affected		
PS 3.2 – Section A.1		
Correction Wording:		

**Table A.1-2
 UID VALUES**

UID Value	UID NAME	Category
...		
1.2.840.10008.1.20.1	Storage Commitment Push Model SOP Class	Workflow Management
<u>1.2.840.10008.1.40</u>	<u>Procedural Event Logging SOP Class</u>	<u>Workflow Management</u>
1.2.840.10008.3.1.2.1.1	Detached Patient Management SOP Class	Workflow Management
...		
1.2.840.10008.3.1.2.6.1	Detached Interpretation Management SOP Class	Workflow Management
<u>1.2.840.10008.4.2</u>	<u>Storage Service Class</u>	<u>Transfer</u>
1.2.840.10008.5.1.1.1	Basic Film Session SOP Class	Print Management
...		
1.2.840.10008.5.1.1.16.376	Printer Configuration Retrieval SOP Class	Print Management

<u>1.2.840.10008.5.1.1.17.376</u>	Printer Configuration Retrieval SOP Instance	Print Management
1.2.840.10008.5.1.1.18	Basic Color Print Management Meta SOP Class	Print Management
...		
1.2.840.10008.5.1.1.32	Pull Stored Print Management Meta SOP Class	Print Management
<u>1.2.840.10008.5.1.1.33</u>	Media Creation Management SOP Class UID	Print Management
1.2.840.10008.5.1.4.1.1.1	Computed Radiography Image Storage	Transfer
...		
1.2.840.10008.5.1.4.1.1.2	CT Image Storage	Transfer
<u>1.2.840.10008.5.1.4.1.1.2.1</u>	Enhanced CT Image Storage	Transfer
1.2.840.10008.5.1.4.1.1.3.1	Ultrasound Multi-frame Image Storage	Transfer
...		
1.2.840.10008.5.1.4.1.1.66	Raw Data Storage	Transfer
<u>1.2.840.10008.5.1.4.1.1.66.1</u>	Spatial Registration Storage	Transfer
<u>1.2.840.10008.5.1.4.1.1.66.2</u>	Spatial Fiducials Storage	Transfer
1.2.840.10008.5.1.4.1.1.77.1.1	VL Endoscopic Image Storage	Transfer
<u>1.2.840.10008.5.1.4.1.1.77.1.1.1</u>	Video Endoscopic Image Storage	Transfer
1.2.840.10008.5.1.4.1.1.77.1.2	VL Microscopic Image Storage	Transfer
<u>1.2.840.10008.5.1.4.1.1.77.1.2.1</u>	Video Microscopic Image Storage	Transfer
1.2.840.10008.5.1.4.1.1.77.1.3	VL Slide-Coordinates Microscopic Image Storage	Transfer
1.2.840.10008.5.1.4.1.1.77.1.4	VL Photographic Image Storage	Transfer
<u>1.2.840.10008.5.1.4.1.1.77.1.4.1</u>	Video Photographic Image Storage	Transfer
<u>1.2.840.10008.5.1.4.1.1.77.1.5.1</u>	Ophthalmic Photography 8 Bit Image Storage	Transfer
<u>1.2.840.10008.5.1.4.1.1.77.1.5.2</u>	Ophthalmic Photography 16 Bit Image Storage	Transfer
<u>1.2.840.10008.5.1.4.1.1.77.1.5.3</u>	Stereometric Relationship Storage	Transfer

1.2.840.10008.5.1.4.1.1.88.11	Basic Text SR	Transfer
1.2.840.10008.5.1.4.1.1.88.22	Enhanced SR	Transfer
1.2.840.10008.5.1.4.1.1.88.33	Comprehensive SR	Transfer
<u>1.2.840.10008.5.1.4.1.1.88.40</u>	<u>Procedure Log Storage</u>	<u>Transfer</u>
1.2.840.10008.5.1.4.1.1.88.50	Mammography CAD SR	Transfer
...		
1.2.840.10008.5.1.4.32	General Purpose Worklist Management Meta SOP Class	Workflow Management
<u>1.2.840.10008.5.1.4.33</u>	<u>Instance Availability Notification SOP Class</u>	<u>Workflow Management</u>
<u>1.2.840.10008.5.1.4.37.1</u>	<u>General Relevant Patient Information Query</u>	<u>Query/Retrieve</u>
<u>1.2.840.10008.5.1.4.37.2</u>	<u>Breast Imaging Relevant Patient Information Query</u>	<u>Query/Retrieve</u>
<u>1.2.840.10008.5.1.4.37.3</u>	<u>Cardiac Relevant Patient Information Query</u>	<u>Query/Retrieve</u>