

DICOM Correction Item

Correction Number		CP-535
Log Summary: Restore mammography view code sequence		
Type of Modification	Name of Standard	
Correction	PS 3.3 2004	
Rationale for Correction		
Whilst editing in the general anatomy macros, the view code sequence and modifier were inadvertently and incorrectly removed from the Mammography Image Module. Also fix the incorrect macro table reference.		
Sections of documents affected		
PS 3.3 C.8.11.7		
Correction Wording:		

Amend PS 3.3 C.8.11.7:

C.8.11.7 Mammography Image Module

...

**Table C.8-74
MAMMOGRAPHY IMAGE MODULE ATTRIBUTES**

Attribute Name	Tag	Type	Attribute Description
...
<i>Include 'General Anatomy Mandatory Macro' Table 10-x-6 10-5</i>			Defined Context ID for the Anatomic Region Sequence is 4013.
<u>View Code Sequence</u>	<u>(0054,0220)</u>	<u>1</u>	<u>Sequence that describes the projection of the anatomic region of interest on the image receptor.</u> <u>Only a single Item shall be permitted in this Sequence.</u>
<u>>Include 'Code Sequence Macro' Table 8.8-1.</u>		<u>Enumerated Value for Context ID is 4014.</u>	
<u>>View Modifier Code Sequence</u>	<u>(0054,0222)</u>	<u>2</u>	<u>View modifier.</u> <u>Zero or more Items may be included in this Sequence.</u>
<u>>>Include 'Code Sequence Macro' Table 8.8-1.</u>		<u>Enumerated Value for Context ID is 4015.</u>	