

### DICOM Correction Item

Correction Number		CP-662
Log Summary: Correction for Enhanced XA/XRF Image Module		
Type of Modification	Name of Standard	
Corrections	PS 3.3 – 2004	
Rationale for Correction		
Correction for Enhanced XA/XRF Image Module:		
Nesting level Patient Gantry Relationship Code Sequence include wrong.		
Referenced Other Plane Sequence includes wrong Macro Table reference.		
Sections of documents affected		
PS 3.3 C.8.19.2		
Correction Wording:		

#### C.8.19.2 Enhanced XA/XRF Image Module

...

**Table C.8.19.2-1**  
**Enhanced XA/XRF Image Module Table**

Attribute Name	Tag	Type	Attribute Description
...			
Patient Gantry Relationship Code Sequence	(0054,0414)	2C	Sequence that describes the orientation of the patient with respect to the head of the table. See Section C.8.4.6.1.3 for further explanation.  Only a single Item shall be permitted in this Sequence.  Required if Positioner Type (0018,1508) equals CARM and C-arm Positioner Tabletop Relationship (0018,9474) equals YES. May be present otherwise.
➤➤Include 'Code Sequence Macro' Table 8.8-1.			Baseline Context ID 21
...			
Referenced Other Plane Sequence	(0008,9410)	1C	A sequence that identifies the SOP Class/Instance pairs of the corresponding plane for a Biplane acquisition device.  Only a single Item shall be permitted in this Sequence.  Required if Image Type (0008,0008) Value 3 is BIPLANE A or BIPLANE B.
➤Include 'SOP Instance Reference Macro' Table C.10-3			
...			

