

## DICOM Correction Item

|  |                  |        |
|--|------------------|--------|
| Correction Number  |                  | CP 702 |
| Log Summary: Author Observer Sequence multiplicity   |                  |        |
| Type of Modification   | Name of Standard |        |
| Clarification  | PS 3.3 2007      |        |
| Rationale for Correction   |                  |        |
| <p>Author Observer Sequence should allow multiple Items to harmonize with HL7 CDA. Also, this and related sequences are Type 3 attributes but are specified to contain zero or more items; normally Type 3 sequences are not defined as being allowed to be empty (the entire attribute is omitted instead of sending an empty attribute).</p> |                  |        |
| Sections of documents affected   |                  |        |
| PS 3.3 C.17.2  |                  |        |
| Correction Wording:  |                  |        |

### C.17.2 SR Document General Module

...

**Table C.17-2**  
**SR DOCUMENT GENERAL MODULE ATTRIBUTES**

| Attribute Name   | Tag         | Type | Attribute Description   |
|--|-------------|------|---|
| ...  | ...         | ...  | ...   |
| Author Observer Sequence   | (0040,A078) | 3    | <p>The person or device that created the clinical content of this document. This attribute sets the default Observer Context for the root of the content tree.</p> <p><del>Zero or one</del> <b>One or more</b> Items may be included in this sequence.</p>   |
| <i>&gt;Include 'Identified Person or Device Macro' Table C.17-3b</i> |             |      |   |
| Participant Sequence   | (0040,A07A) | 3    | <p>Persons or devices related to the clinical content of this document.</p> <p><del>Zero</del> <b>One</b> or more Items may be included in this sequence.</p>   |
| >Participation Type  | (0040,A080) | 1    | <p>Participant's role with respect to the clinical content of this document. See C.17.2.5.</p> <p>Defined Terms:</p> <p>SOURCE – Equipment that contributed to the content</p> <p>ENT – Data enterer (e.g., transcriptionist)</p> <p>ATTEST – Attestor</p> <p>Note: In HL7 Structured Documents, the participation comparable to Attestor is the "Authenticator".</p> |
| >Participation Datetime  | (0040,A082) | 2    | Datetime of participation with respect to the   |

|  |             |     |   |
|--|-------------|-----|---|
|  |             |     | clinical content of this document.  |
| <i>&gt;Include 'Identified Person or Device Macro' Table C.17-3b</i> |             |     |   |
| Custodial Organization Sequence                                      | (0040,A07C) | 3   | <p>Custodial organization for this SR Document instance. Represents the organization from which the document originates and that is in charge of maintaining the document, i.e., the steward of the original source document.</p> <p>Note: This may or may not be identical to the Institution identified in the Equipment Module.</p> <p><b><u>Zero or one Items may be included Only a single item shall be present</u></b> in this sequence.</p> |
| <i>&gt;Institution Name</i>  | (0008,0080) | 2   | Name of Custodial Institution or Organization.  |
| <i>&gt;Institution Code Sequence</i>                                 | (0008,0082) | 2   | <p>Coded identifier of Custodial Institution or Organization.</p> <p>Zero or one Items may be included in this sequence.</p>  |
| <i>&gt;&gt;Include 'Code Sequence Macro' Table 8.8-1</i>             |             |     | <i>No Baseline Context ID defined.</i>  |
| ...  | ...         | ... | ...   |