

DICOM Correction Item

Correction Number		CP 703
Log Summary: Include Series Description, Date, Time in SR IODs		
Type of Modification	Name of Standard	
Addition	PS 3.3 2007	
Rationale for Correction		
The SR-related IODs do not include the General Series Module, and do not specify the presence of Series Description, Series Date and Series Time, which is very useful for all browsers.		
Sections of documents affected		
PS 3.3 C.17.1, C.17.6.1		
Correction Wording:		

C.17.1 SR Document Series Module

Table C.17-1 defines the Attributes of the SR Document Series. A Series of SR Documents may contain any number of SR Documents.

Note: Series of SR Documents are separate from Series of Images or other Composite SOP Instances. SR Documents do not reside in a Series of Images or other Composite SOP Instances.

Table C.17-1
SR DOCUMENT SERIES MODULE ATTRIBUTES

Attribute Name	Tag	Type	Attribute Description
Modality	(0008,0060)	1	Modality type. Enumerated Value: SR = SR Document
Series Instance UID	(0020,000E)	1	Unique identifier of the Series. Note: No SR-specific semantics are specified.
Series Number	(0020,0011)	1	A number that identifies the Series. Note: No SR-specific semantics are specified.
<u>Series Date</u>	<u>(0008,0021)</u>	<u>3</u>	<u>Date the Series started.</u>
<u>Series Time</u>	<u>(0008,0031)</u>	<u>3</u>	<u>Time the Series started.</u>
<u>Series Description</u>	<u>(0008,103E)</u>	<u>3</u>	<u>User provided description of the Series</u>
Referenced Performed Procedure Step Sequence	(0008,1111)	2	Uniquely identifies the Performed Procedure Step SOP Instance for which the Series is created. Only a single Item shall be permitted in this sequence. Notes: 1. The Performed Procedure Step referred to by this Attribute is the Step during which this Document is generated. 2. If this Document is generated during the same Modality or General Purpose Performed Procedure Step as the evidence in the current interpretation procedure, this attribute may

			<p>contain reference to that Modality or General Purpose Performed Procedure Step.</p> <p>3. This Attribute is not used to convey reference to the evidence in the current interpretation procedure. See Current Requested Procedure Evidence Sequence (0040,A375).</p> <p>4. This Sequence may be sent zero length if the Performed Procedure Step is unknown.</p>
>Referenced SOP Class UID	(0008,1150)	1C	<p>Uniquely identifies the referenced SOP Class.</p> <p>Required if Referenced Performed Procedure Step Sequence (0008,1111) is sent.</p>
> Referenced SOP Instance UID	(0008,1155)	1C	<p>Uniquely identifies the referenced SOP Instance.</p> <p>Required if Referenced Performed Procedure Step Sequence (0008,1111) is sent.</p>

C.17.6.1 Key Object Document Series Module

Table C.17.6-1 defines the Attributes of the Key Object Document Series.

Note: Series of Key Object Selection Documents are separate from Series of Images or other Composite SOP Instances. Key Object Documents do not reside in a Series of Images or other Composite SOP Instances.

**Table C.17.6-1
 KEY OBJECT DOCUMENT SERIES MODULE ATTRIBUTES**

Attribute Name	Tag	Type	Attribute Description
Modality	(0008,0060)	1	<p>Modality type.</p> <p>Enumerated Value:</p> <p>KO = Key Object Selection</p>
Series Instance UID	(0020,000E)	1	<p>Unique identifier of the Series.</p> <p>Note: No specific semantics are specified.</p>
Series Number	(0020,0011)	1	<p>A number that identifies the Series.</p> <p>Note: No specific semantics are specified.</p>
<u>Series Date</u>	<u>(0008,0021)</u>	<u>3</u>	<u>Date the Series started.</u>
<u>Series Time</u>	<u>(0008,0031)</u>	<u>3</u>	<u>Time the Series started.</u>
<u>Series Description</u>	<u>(0008,103E)</u>	<u>3</u>	<u>User provided description of the Series</u>
Referenced Performed Procedure Step Sequence	(0008,1111)	2	<p>Uniquely identifies the Performed Procedure Step SOP Instance for which the Series is created. Only a single Item shall be permitted in this sequence.</p> <p>Notes: See notes on this attribute in Section C.17.1 SR Document Series Module</p>
>Referenced SOP Class UID	(0008,1150)	1C	<p>Uniquely identifies the referenced SOP Class.</p> <p>Required if Referenced Performed Procedure Step Sequence (0008,1111) is sent.</p>
> Referenced SOP Instance UID	(0008,1155)	1C	<p>Uniquely identifies the referenced SOP Instance.</p> <p>Required if Referenced Performed Procedure Step Sequence (0008,1111) is sent.</p>

