

DICOM Correction Item

Correction Number		CP-862
Log Summary: Ophthalmic Frame Location Functional Group Macro is unused		
Type of Modification	Name of Standard	
Correction	PS 3.3 2008	
Rationale for Correction		
The Ophthalmic Tomography IOD does not include the Ophthalmic Frame Location Functional Group Macro in the table of Functional Group Macros.		
Sections of documents affected		
PS 3.3 A.52.4.3		
Correction Wording:		

Amend PS 3.3 A.52.4.3:

A.52.4.3 Ophthalmic Tomography Image Functional Group Macros

Table A.52.4.3-1 specifies the use of the Functional Group macros used in the Multi-frame Functional Group Module for the Ophthalmic Tomography Image IOD.

**Table A.52.4.3-1
 OPHTHALMIC TOMOGRAPHY FUNCTIONAL GROUP MACROS**

Function Group Macro	Section	Usage
...
Contrast/Bolus Usage	C.7.6.16.2.12	C – Required if Contrast/Bolus Agent Sequence (0018,0012) is used. May not be used as a Shared Functional Group
<u>Ophthalmic Frame Location</u>	<u>C.8.17.10.1</u>	<u>U</u>